

Q BANDAGES

QUALITY GERMAN MANUFACTURE

Q HOSPITAL CREPE

- For general purpose use
- Lightweight open fabric composition
- Non fraying edges
- Economical
- Available as single rolls and Packs of 12



Cotton Crepe

Q CONFORMING

- For all types of dressing retention
- Soft & comfortable on the skin
- Conforms well to body contours without creasing
- Available in Packs of 12



Soft & Comfortable

Q MEDIUM CREPE

- Medium weight
- Wrinkled crepe
- Washable cotton
- Non-fraying edges
- Light support
- Available as single rolls and Packs of 12



Light Support

Q COHESIVE RETENTION

- For dressing retention
- Self adhesive
- Doesn't stick to skin
- Easy & time saving application
- Economical
- Lightweight & secure dressing retention for high movement areas
- Available as single rolls and Packs of 12



4m or 20m Lengths

Q HEAVY CREPE

- Moderate compression
- Flesh coloured
- Heavy weight for medium to long term bandaging
- Available as single rolls and Packs of 12



Flesh Colour

Q COHESIVE SUPPORT

- For non-slip support
- Self adhesive & flexible
- Easy pain free removal
- Light support
- Flesh coloured
- Available as single rolls and Packs of 12



Non-Slip

SEE REVERSE FOR ORDERING DETAILS

Q BANDAGES

QUALITY GERMAN MANUFACTURE

ORDERING INFORMATION

RETENTION

RE-ORDER
CODE

Q-Hospital crepe 5cm x 1.5m (unstretched)	Roll	1605
Q-Hospital crepe 7.5cm x 1.5m (unstretched)	Roll	1675
Q-Hospital crepe 10cm x 1.5m (unstretched)	Roll	1610
Q-Hospital crepe 15cm x 1.5m (unstretched)	Roll	1615
Q-Medium crepe 2.5cm x 1.6m (unstretched)	Roll	1502
Q-Medium crepe 5cm x 1.6m (unstretched)	Roll	1505
Q-Medium crepe 7.5cm x 1.6m (unstretched)	Roll	1507
Q-Medium crepe 10cm x 1.6m (unstretched)	Roll	1510
Q-Medium crepe 15cm x 1.6m (unstretched)	Roll	1515
Q-Conforming elastic retention 2.5cm x 1.5m (unstretched)	Pack 12	1225
Q-Conforming elastic retention 5cm x 1.5m (unstretched)	Pack 12	1250
Q-Conforming elastic retention 7.5cm x 1.5m (unstretched)	Pack 12	1275
Q-Conforming elastic retention 10cm x 1.5m (unstretched)	Pack 12	1210
Q-Conforming elastic retention 15cm x 1.5m (unstretched)	Pack 12	1215
Q-Cohesive elastic retention 2.5cm x 4m (stretched)	Roll	1902
Q-Cohesive elastic retention 6cm x 4m (stretched)	Roll	1906
Q-Cohesive elastic retention 6cm x 20m (stretched)	Roll	1907
Q-Cohesive elastic retention 8cm x 4m (stretched)	Roll	1908
Q-Cohesive elastic retention 8cm x 20m (stretched)	Roll	1909
Q-Cohesive elastic retention 10cm x 4m (stretched)	Roll	1910
Q-Cohesive elastic retention 10cm x 20m (stretched)	Roll	1911
Q-Cohesive elastic retention 12cm x 4m (stretched)	Roll	1912

SUPPORT

Q-Heavy crepe 5cm x 2.3m (unstretched)	Roll	1405
Q-Heavy crepe 7.5cm x 2.3m (unstretched)	Roll	1475
Q-Heavy crepe 10cm x 2.3m (unstretched)	Roll	1410
Q-Heavy crepe 15cm x 2.3m (unstretched)	Roll	1415
Q-Cohesive Support elastic non woven 2.5cm x 4.5m (stretched)	Roll	1802
Q-Cohesive Support elastic non woven 5cm x 4.5m (stretched)	Roll	1825
Q-Cohesive Support elastic non woven 7.5cm x 4.5m (stretched)	Roll	1806
Q-Cohesive Support elastic non woven 10cm x 4.5m (stretched)	Roll	1808
Q-Cohesive Support elastic non woven 15cm x 4.5m (stretched)	Roll	1810



SURGI SUPPLIES PTY LTD

3/29 MACQUARIE DRIVE
THOMASTOWN VIC. 3074

TEL: (03) 9466 2803
(03) 9464 5691

FAX: (03) 9464 0237

EMAIL: sales@surgisupplies.com.au

Australasian Agent: Healthmed Marketing Australia • Ph: (03) 9466 2697