



Designed to optimize the efficiency and quality of dental care in healthcare facilities

- Robustness: excellent resistance, natural anodised structure, ABS-mass tinted side walls
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ergonomic and handy: designed for medical staff comfort



Robust



Hygienic



Ergonomic



Customizable

Colour

White

- Single-column trolley
- Height: 100 cm
- Roller shutter door
- Code lock (with spare key)
- 1 push handle
- 3 standard rails for accessories
- Delivered without equipment

100 cm height trolley*

8Y144.DT

This trolley can be equipped with five drawers for the storage of dental prostheses cleaning boxes.



Accessories and internal fittings

Extension shelf, standard rail, arm for laptop, upper frame, waste bin, drawer...

Please refer to product data sheets (please use compatible dental care trolley accessories)



Technical data



* Usable dimensions

** Usable height on work surface

Materials

Structure	Anodised aluminium profile
Work surface	ABS / PMMA thermoplastic
Slides	ABS (Acrylonitrile Butadiene Styrene)
Closing system	Roller shutter door: grey PP (Polypropylene)
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)

Additional Information

Wheels	100 mm diameter– 4 swivel wheels including 1 with brake - wheel and wheel-pin double blocking
Empty weight	21 kg
Loading capacity	5 kg – 70 mm levels
Capacity	10 levels of 70 mm
Locking system	Code lock (with spare key)
External side panel colour	ABS (Acrylonitrile Butadiene Styrene) white
Equipment	2 standard rails on the right side 1 rear standard rail 1 push handle

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions L X D x H	80 x 120 x 120 cm
Gross weight (trolley + packaging)	31 kg	Volume	1 m ³