

EZbra[®] Breast Dressing

About EZbra[®]

EZbra[®] represents an elegant solution specifically crafted to cater to all the post-breast surgery needs of women. It stands as a sterile, single-use, medical-grade compression bra, while simultaneously serving as a non-adherent and absorbent primary wound dressing solution.

EZbra[®] features three distinct and individually adjustable layers of compression, ensuring optimal comfort and support during the recovery process.


1 out of 8 women will be diagnosed with breast cancer during her lifetime



About **17,000,000** breast procedures are performed annually worldwide.



**Cancer
Surgeries**



**Cosmetic
Surgeries**



Biopsies

And yet, there is currently no established standard of care for post-operative breast dressings.

Everyone Benefits From EZbra[®]

Patients

Single accessory; Intuitive dressing; Easy-to-change ; Comfortable ; Independent use ; Physical / Psychological support; Fast recovery; Maximum freedom & independence;

Surgeons

Fully-integrated sterile post-op dressing; Implant stabilization; Compression; Absorption; Individually adjustable breast control; Support for drain management;

Nurses

No-leak; Short learning curve; Fast & easy replacement; Can be applied by patients; Great patient compliance and reviews;

>> EZbra is CE-Marked, FDA Class I-designated (exempt from 510K), ISO certified.

EZbra[®] Breast Dressing

All-In-One Single-Use
Wearable Solution

>> Contact Us <<

Web: ezbra.net

Mail: info@ezbra.net



**Controlled
Compression**

Provides needed
compression & support.
Hygienic solution that
won't compromise
support needs.



Absorption

Internal absorbent
layer eliminates the
need for absorbing
pads, gauze & tapes.



**Complies
With Asymmetry**



**Drain
Management**

Drain tube
& bulb holder.



**Implant
Stabilization**

Integrated upper & lower
compression bands.



**Easy
Application**



**Learn More
About EZbra[®]**



What is
EZbra[®] Breast Dressing?



Patients
Instructions for use



Medical Team
Instructions for use

**Select Your
Perfect Fit**

Our simple size chart can help you find
the right EZbra for maximum comfort

US	28	30	32	34	36	38	40	42	44	46	48
EU	60	65	70	75	80	85	90	95	100	105	110
A	XS	XS	XS	S	S	S	M	M	M	M	L
B	XS	XS	XS	S	S	S	M	M	M	M	L
C	XS	XS	S	S	S	M	M	M	M	L	L
D	XS	S	S	S	S	M	M	M	M	L	L
DD	S	M	M	M	M	M	M	L	L		
DDD		M	M	M	M	M	M	L	L		