

# Wondaleaf® UriFYX Urinary Catheter

Transparent and breathable



Wondaleaf UriFYX U-Series and R-Series are adhesive, free-size, and breathable external urinary catheters designed to stay securely on the penis for long periods. The U-Series is suitable for both circumcised and uncircumcised, while the R-Series features a reservoir that can be emptied as needed.

Made from soft, medical-grade polyurethane and free from latex, they prioritize comfort and safety. Their transparent nature allows for easy detection of infections, and the conformable design adjusts to changes in penis size. Compatible with most urine bags, ensure the bag is below the penis height for easy, mess-free urine management.



Free-Size

Fits and conforms to any change in penis size



Breathable

Prevents irritation and moisture buildup



Transparent

Allows visual monitoring of penis

## Stick-on external urinary catheter

U-Series



9555864005070

R-Series



9555864005049

### Information

Wondaleaf UriFYX is suitable for:

- Urinary Incontinence:** Unintentional urine leakage due to loss of bladder control.
- Overactive Bladder:** Urgent urination with difficulty reaching the bathroom in time.
- Limited Mobility:** Ideal for those who are bedridden or have difficulty accessing a bathroom.
- Dementia:** Helps individuals who struggle with recognizing the need to urinate.

Hold up to  
**400cc of urine!**



MDA Reg No: GA7462924-178491



Urinary Problems



Connect to a urine bag



Comfortably packed inside underwear!

## Non-invasive method Condom Catheter

## Current method Foley Ballon Catheter Painful Indwelling Catheter

### A Zero Pain Solution in Male Urine Management

Suitable for Male Patients suffering from  
• Dementia • Urinary incontinence • Bedridden



Fits and conforms to any change in penis size



Prevents irritation and moisture buildup



Allows visual monitoring of penis



Stick-to-use



Risk of Inflammation

