

Moulding The Future of Medical Plasti



We work with a wide range of industry sectors including:



#### STERILE SERVICES & DECONTAMINATION

Operating Room,  
Instrument Packing



#### PERSONAL USE & HOME CARE

Daily Living Aids,  
Community Nurse,  
Incontinent Care



#### HOSPITALS

Surgery, Clinics, Wards,  
Veterinary, Dentistry



#### NURSING HOMES

Patient Care, Rehabilitation



#### CUSTOM PACK MANUFACTURERS

Single Use Packs, Contract  
Instrument Processing

### WARWICK - reusable range

Compartment Trays	page	6
Disinfectant & Soaking Containers	page	6
Instrument Trays	page	7-8
DIN Style Baskets	page	8
Instrument Tray Lids	page	9
Gallipots	page	9
Lotion & Sponge Bowls	page	10
Jugs	page	10
Kidney Dishes	page	11
Theatre Bowls	page	11
Wash Bowls	page	11
Medical Boxes, Instrument Boxes & Storage Containers	page	12
Scope Containers	page	13
Silicone Protection	page	13
Tray Tags	page	14

### STERILITE - single use range

Kidney Dishes	page	16
Compartment Trays	page	17
Gallipots	page	17
Quivers	page	18
Sponge Bowls	page	18
Theatre Bowls	page	18

### SPA - care and living aid range

Bedpans	page	20-21
Bedpan Holders	page	21
Commode Pans	page	21
Bidets	page	21
Urinal Bottles and Urinal Pans	page	22
Urinal Accessories	page	22
Drinking Beakers and Cups	page	23
Drinking Cup Lids	page	23
Denture Cups	page	23
Jugs and Jug Sets	page	24
Medicine Measures	page	24
Patient Trays	page	25
Kidney Dishes	page	25
Vomit Bowls	page	25



To find out more about our service, product ranges, product decontamination cycle and industry news:

visit [www.sasco.co.uk](http://www.sasco.co.uk)



# WARWICK RANGE

## Reusable Medical Polyware

The Warwick Range is a wide range of sterilizable long-life plastic polyware products.

Designed as an alternative to stainless steel products, they are ideal for use in Decontamination Services, Operating Rooms, Clinics, Surgeries and other surgical or healthcare services.

Products are designed for repeated washing, disinfection and sterilization. The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>). Products comply with (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.

## Warwick Benefits

- Low cost.
- Robust and hard wearing - guaranteed for 1000 cycles.
- Easy to clean - blood protein is easily removed from the surface, leaving no residue.
- 50% lighter than stainless steel, assisting users with storage, lifting and handling.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>).
- Can be used as standalone items or as a support component in surgical sets.
- Widely used in Sterile Services for packing instrument sets and procedure kits for the autoclave.
- Disinfection sets can be used for manual sterilization in Sterile Services, Clinics, etc.
- Used non-sterile in hospital wards and clinics.
- Smooth corners and surfaces, reducing risk of paper rips and the need for re-processing.
- Autoclave sterilization temperature is 134.6°C + 2°C (generally referred to as 135°C).



# Compartment Trays, Disinfectant and Soaking Containers



Compartment Tray -  
3 compartments

**Code: CIT1814**

Dimensions: 181x134x54mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 38



Compartment Tray -  
3 compartments

**Code: CIT2013**

Dimensions: 205x130x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 120



Compartment Tray -  
4 compartments

**CIT2718**

Dimensions: 270x180x41mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 25



Compartment Tray -  
4 compartments

**Code: LCT2215**

Dimensions: 220x150x24mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 160



Disinfection Set -  
Soaking container, strainer  
and hinged lid  
**Code: DSSW1000**

Dimensions: 260x130x65mm  
Volume: 1 litre  
Material: Polypropylene  
Colour: White  
Case Size: 10



Disinfection Set -  
Soaking container, strainer  
and lid  
**Code: DSS2711**

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: Blue  
Case Size: 5



Disinfection Set -  
Soaking container, strainer  
and lid  
**Code: DSSW2711**

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: White  
Case Size: 5



Disinfection Set -  
Soaking container, strainer  
and vented lid  
**Code: DSSV2711**

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: Blue  
Case Size: 5



Disinfection Set -  
Soaking container, strainer,  
drip tray and lid  
**Code: DSS4215**

Dimensions: 426x156x138mm  
Volume: 4litre  
Material: Polypropylene  
Colour: Blue  
Case Size: 1



Disinfection Set -  
Soaking container, strainer,  
drip tray and lid  
**Code: DSS8515**

Dimensions: 890x200x145mm  
Volume: 12litre  
Material: Polypropylene  
Colour: Blue  
Case Size: 1



Disinfection Set - 2 x soaking con-  
tainers, 1 x strainer, 2 x lids and 1  
x drip tray  
**Code: DSS8515TP**

Dimensions: 890x200x145mm  
Volume: 12litre  
Material: Polypropylene  
Colour: Blue  
Case Size: 1



**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free



# Instrument Trays and DIN Style Baskets



Instrument Tray -  
Solid base

**Code: IT1813**

Dimensions: 181x134x54mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 40

R W S E HP LF



Instrument Tray -  
Solid base

**Code: IT1913**

Dimensions: 195x135x23mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 115

R W S E HP LF



Instrument Tray -  
Solid base

**Code: T2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 60

R W S E HP LF



Instrument Tray -  
Solid ribbed base

**Code: IT2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 60

R W S E HP LF



Instrument Tray -  
Solid base

**Code: T3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 25

R W S E HP LF



Instrument Tray -  
Solid base

**Code: IT2715**

Dimensions: 270x150x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 60

R W S E HP LF

Distributed By



Instrument Tray -  
Solid ribbed base

**Code: IT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 25

R W S E HP LF



Instrument Tray -  
Solid ribbed base

**Code: IT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 12

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

# Instrument Trays and DIN Style Baskets



Instrument Tray -  
Mesh base and sides

**Code: MIT1515**

Dimensions: 200x180x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 50

R W S E HP LF



Instrument Tray -  
Mesh base

**Code: MIT2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 60

R W S E HP LF



Instrument Tray -  
Mesh base and sides

**Code: MIT2709**

Dimensions: 275x110x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 60

R W S E HP LF



Instrument Tray -  
Mesh base

**Code: MIT2715**

Dimensions: 270x150x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 60

R W S E HP LF



Instrument Tray -  
Mesh base and sides

**Code: MIT2817**

Dimensions: 280x170x55mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 25

R W S E HP LF



Instrument Tray -  
Mesh base and sides

**Code: MIT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 24

R W S E HP LF



Instrument Tray -  
Mesh base

**Code: MIT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 12

R W S E HP LF



Instrument Tray -  
Perforated base

**Code: PIT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 25

R W S E HP LF



Instrument Tray -  
Perforated base

**Code: PIT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 12

R W S E HP LF



DIN Style Tray -  
Mesh base and sides

**Code: DINT2525**

Dimensions: 255x255x50mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 20

R W S E HP LF



DIN Style Tray -  
Mesh base and sides

**Code: DINT4825**

Dimensions: 480x255x50mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 25

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization


E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free



# Instrument Tray Lids



Instrument Tray Lid -  
For trays IT1813/CIT1814

**Code: ITL181**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



Instrument Tray Lid -  
For trays IT2015/MIT2015/T2015  
**Code: ITL201**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



Instrument Tray Lid -  
For trays IT2715/MIT2715

**Code: ITL271**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



Instrument Tray Lid -  
For tray CIT2718

**Code: ITL281**


Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



Instrument Tray Lid -  
For tray MIT2817

**Code: ITL287**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



Instrument Tray Lid -  
For trays IT3025/MIT3025/  
PIT3025/T3025  
**Code: ITL302**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



# Gallipots



Gallipot

**Code: GP40**

Dimensions: 40x30mm  
Volume: 35ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 100



Gallipot

**Code: GP60**

Dimensions: 60x40mm  
Volume: 95ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 300



Gallipot

**Code: GP80**

Dimensions: 80x60mm  
Volume: 280ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 80



lipot - Flat base  
iduated

**Je: GPLB60**

Dimensions: 60x40mm  
Volume: 50ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 400



Gallipot - Flat base  
Graduated

**Code: GPLB80**

Dimensions: 80x60mm  
Volume: 200ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 200



# Lotion and Sponge Bowls

Distributed By



Lotion Bowl - Flat base  
Graduated

Code: GPLB100

Dimensions: 100x45mm  
Volume: 300ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 135



Lotion Bowl - Flat base  
Graduated

Code: GPLB150

Dimensions: 150x70mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 60



Sponge Bowl - Flat base  
Graduated

Code: SB140

Dimensions: 140x65mm  
Volume: 300ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 160



Lotion Bowl  
Graduated

Code: LB100

Dimensions: 100x45mm  
Volume: 280ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 200



Lotion Bowl  
Graduated

Code: LB150

Dimensions: 150x70mm  
Volume: 900ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 60



Lotion Bowl  
Graduated

Code: LB200

Dimensions: 200x90mm  
Volume: 2000ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 30



Lotion Bowl  
Graduated

Code: LB250

Dimensions: 250x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 20



# Jugs



Jug  
Graduated

Code: JG500

Dimensions: 95mm x 95mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 50



Jug  
Graduated

Code: JG1000

Dimensions: 120mm x 125mm  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 45



Jug  
Graduated

Code: JG2000

Dimensions: 130mm x 200mm  
Volume: 2000ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 20





# Kidney Dishes, Theatre Bowls and Wash Bowls



Kidney Dish  
Graduated

**Code: KD150**

Dimensions: 150x35mm  
Volume: 150ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 210



Kidney Dish  
Graduated

**Code: KD200**

Dimensions: 200x45mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 150



Kidney Dish  
Graduated

**Code: KD250**

Dimensions: 250x55mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 60



Kidney Dish  
Graduated

**Code: KDY250**

Dimensions: 250x55mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Yellow  
Case Size: 60



Kidney Dish  
Graduated

**Code: KD300**

Dimensions: 300x65mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 25



Kidney Dish Tray -  
Shallow kidney shaped tray

**Code: KDT185**

Dimensions: 185x92mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 70



Wash Bowl  
Graduated

**Code: WBN3112**

Dimensions: 310x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 25



Wash Bowl  
Graduated

**Code: WBN3513**

Dimensions: 350x120mm  
Volume: 5000ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 25



Theatre Bowl  
Graduated

**Code: WB3112**

Dimensions: 310x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 25



Theatre Bowl  
Graduated

**Code: WB3513**

Dimensions: 350x120mm  
Volume: 5000ml  
Material: Polypropylene  
Colour: Blue  
Case Size: 25



**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free

# Medical Boxes, Instrument Boxes & Storage Containers



Instrument Box -  
Storage box with lid

Code: IBX2718

Dimensions: 257x180x128mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



Instrument box -  
Storage box with vented lid

Code: IBXV2718

Dimensions: 257x180x128mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 10



Instrument Box -  
Storage box with lid

Code: IBXN290

Dimensions: 300x224x140mm  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Case Size: 12



Distributed By



Transportation Box -  
With hinged lid

Code: MB2318

Dimensions: 230x170x200mm  
Volume: N/A  
Material: ABS  
Colour: Orange  
Other colours available  
Case Size: 10



- Reusable
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free



## Scope Containers



Scope Container -  
With silicone dividers and lid

**Code: SMT5510**

Dimensions: 660x110x66mm  
Volume: N/A  
Material: Polypropylene  
Colour: Lilac  
Case Size: 1



Silicone Plug -  
For SMT5510

**Code: SG5510**

Dimensions: N/A  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 10




## Silicone Protection



Silicone Pad -  
For trays IT/MIT2715 and  
small metal DIN tray  
**Code: SP2313**

Dimensions: 230x132mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 10




Silicone Pad -  
For trays IT/MIT/PIT/T3025  
**Code: SP2823**

Dimensions: 280x230mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 10



Silicone Pad -  
Fits 1/2 tray IT/MIT/PIT4030  
**Code: SP2719**

Dimensions: 270x190mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 10



Silicone Pad -  
Fits 1/2 tray DINT4825  
**Code: SP4622**

Dimensions: 230x220mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 10




Silicone Pad -  
For medium metal DIN tray  
**Code: SP2324**

Dimensions: 235x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 10



Silicone Pad -  
For large metal DIN tray  
**Code: SP4724**

Dimensions: 470x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 10



Silicone Roll -  
1 metre length with 10  
divisible sections  
**Code: SR10**

Dimensions: 100x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Case Size: 1



# Tray Tags



Tray Tag

**Code: TAG-010**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: White  
Case Size: 100



Tray Tag

**Code: TAG-015**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Yellow  
Case Size: 100



Tray Tag

**Code: TAG-020**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Blue  
Case Size: 100



Tray Tag

**Code: TAG-022**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Light Green  
Case Size: 100



Tray Tag

**Code: TAG-025**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Orange  
Case Size: 100



Tray Tag

**Code: TAG-030**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Red  
Case Size: 100



Tray Tag

**Code: TAG-050**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Green  
Case Size: 100



Tray Tag

**Code: TAG-055**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Grey  
Case Size: 100



Tray Tag

**Code: TAG-060**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Lavender  
Case Size: 100



Tray Tag

**Code: TAG-065**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Tan  
Case Size: 100



Tray Tag

**Code: TAG-070**

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Black  
Case Size: 100



Stainless Steel Attachment Ring

**Code: TR30**

Dimensions: 30mm  
Volume: N/A  
Material: Steel  
Colour: Silver  
Case Size: 100



- R Reusable
- W Washer Decontamination
- S Suitable for Steam Sterilization
- E Suitable for Ethylene Oxide Sterilization
- HP Suitable for Hydrogen Peroxide Sterilization
- LF Latex Free

# STERILITE RANGE

## Single Use Medical Polyware

The Sterilite Range is a wide range of single use plastic polyware products.

Designed to speed up packing and sterilizing procedures. Sterilite is ideal for use in Sterile Pack Manufacturing, Decontamination Services, Operating Rooms and Clinics due to its low cost, product design and performance in hospitals.

The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>).

Products comply with (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable. Material used is polypropylene unless otherwise stated.

## Sterilite Benefits

- Low cost and disposable.
- No need to invest time and money into retrieval and decontamination after use.
- Clean and ready for sterilization.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>).
- Can be used as stand-alone items or as a support (supplementary) component in operating rooms and surgical procedures.
- Widely used in Sterile Services for packing instrument sets and procedure kits.
- Used non-sterile in hospital wards and clinics.





# Kidney Dishes



Kidney Dish  
Graduated

**Code: KB250**

Dimensions: 252x100x48mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 450



Kidney Dish  
Graduated

**Code: KB300**

Dimensions: 300x140x57mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 135



Single Use

Washer Decontamination

Suitable for Steam Sterilization

Suitable for Ethylene Oxide Sterilization

Suitable for Hydrogen Peroxide Sterilization

Latex Free

# Compartment Tray



Packing Tray -  
4 compartments

**Code: LCT2215S**

Dimensions: 220x150x24mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Case Size: 160



# Gallipots



Gallipot  
Graduated

**Code: GPS60**

Dimensions: 62x30mm  
Volume: 60ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 2000



Gallipot  
Graduated

**Code: GPS120**

Dimensions: 85x30mm  
Volume: 120ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 1000



Gallipot  
Graduated

**Code: GP100**

Dimensions: 100x53mm  
Volume: 250ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 1050



Gallipot  
Graduated

**Code: GP130**

Dimensions: 130x66mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 735



Single Use

Washer Decontamination

Suitable for Steam Sterilization

Suitable for Ethylene Oxide Sterilization

Suitable for Hydrogen Peroxide Sterilization

Latex Free

# Quiver



**Quiver**

**Code: QV30**

Dimensions: 190x65mm  
 Volume: N/A  
 Material: Polypropylene  
 Colour: Natural  
 Case Size: 75



**Quiver Extension –  
For QV30**

**Code: QV47**

Dimensions: 75x244mm  
 12mm(loop)  
 Volume: N/A  
 Material: Polypropylene  
 Colour: Pink  
 Case Size: 60



# Sponge Bowls and Theatre Bowl

**Sponge Bowl  
Graduated**

**Code: SB170**

Dimensions: 170x84mm  
 Volume: 1000ml  
 Material: Polypropylene  
 Colour: Natural  
 Case Size: 340



**Sponge Bowl  
Graduated**

**Code: SB200**

Dimensions: 210x90mm  
 Volume: 2000ml  
 Material: Polypropylene  
 Colour: Natural  
 Case Size: 220



**Theatre Bowl  
Graduated**

**Code: TB350**

Dimensions: 352x116mm  
 Volume: 5000ml  
 Material: Polypropylene  
 Colour: Blue  
 Case Size: 50



Distributed By

Single Use

Washer Decontamination

Suitable for Steam Sterilization

Suitable for Ethylene Oxide Sterilization

Suitable for Hydrogen Peroxide Sterilization

Latex Free



# SPA RANGE

## Nursing & Care, Daily Living Aids and Accessories

The Spa range includes reusable products for personal living, patient care, toileting and bedside comfort.

The Spa Range is designed for use in Hospital Wards, Rehabilitation, Nursing Homes, Community Care and for Personal Use in the Home.

Products comply with BS2588 For Portable Sanitary Products, BS3215 for Reusable Urinals, (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.

## Spa Benefits

- Toileting products that preserve patient dignity and independence by empowering patients to toilet themselves with little or no assistance.
- Urinal Bottles specifically designed to the male and female anatomy for a comfortable and safe fit during use.
- Product design reduces risk of cross infection, unpleasant odours and sights.
- Easy to use, clean and disinfect,
- meaning both patient and carer can be confident that products are safe to re-use.
- Beakers and Drinking Cups are wide, easy to grip, non-scratch, drop-proof and have a nondrip rim.
- Handled beakers are ideal for patients with unsteady or infirm hands.
- Medicine Measures are sturdy, ideal for long-term use and regularly checked by Trading Standards for measurement accuracy and compliance.
- All Spa products can be cleaned in automatic washer disinfectors. Also suitable for hand washing using detergents and disinfectants.



# Bedpans

Distributed By

Bedpan -  
Adult

Code: BPN

Dimensions: 370x310x90mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Case Size: 30

R W S E HP LF



Slipper Pan -  
Bedpan with handle and lid

Code: FPN

Dimensions: 449x309x99mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Case Size: 12

R W S E HP LF



Slipper Pan -  
Bedpan with handle and lid

Code: FPNW

Dimensions: 449x309x99mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: White  
Case Size: 12

R W S E HP LF



Bedpan -  
Hospital style bedpan

Code: VP

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Case Size: 18

R W S E HP LF



Bedpan -  
Hospital bedpan with lid

Code: VPN

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Case Size: 15

R W S E HP LF



Bedpan -  
Hospital bedpan with lid

Code: VPNW

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: White  
Case Size: 15

R W S E HP LF



Bedpan -  
European bedpan with handle  
and lid

Code: DPNF

Dimensions: 420x315x100mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Case Size: 15

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

## Bedpans and Bedpan Holders

Slipper Pan -  
Female urinal and  
infants bedpan  
Code: SPN

Dimensions: 335x245mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Case Size: 25

R W S E HP LF



Bedpan Support -  
For pulp disposable midi  
slipper pan 1.3 litres  
Code: PPS

Dimensions: 470x290x75mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Cream  
Case Size: 30

R W S E HP LF

Bedpan Support -  
For pulp disposable bedpan  
2 litres  
Code: PBPS

Dimensions: 370x310x90mm  
Volume: 2000ml  
Material: Polypropylene  
Colour: Cream  
Case Size: 30

R W S E HP LF



## Commode Pans and Bidets



Commode Pan / Bidet

Code: CMP

Dimensions: 359x338x130mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Green  
Case Size: 20

R W S E HP LF



Commode Pan / Bidet

Code: CMPW

Dimensions: 359x338x130mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: White  
Case Size: 20

R W S E HP LF



Commode Pan -  
Commode pan with lid

Code: CBL

Dimensions: 259x259x144mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: White  
Case Size: 25

R W S E HP LF

Distributed By

2 Single Use

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free



# Urinal Bottles and Urinal Pans

Urinal Bottle -  
Male urinal with handle  
Graduated  
Code: UML1000

Dimensions: 180x100x85mm  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Case Size: 50

R W S E HP LF



Urinal Bottle -  
Male urinal, hospital style  
Graduated  
Code: UM1000

Dimensions: 170x90x 80mm  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 25

R W S E HP LF



Urinal Bottle -  
Male urinal with handle  
Graduated  
Code: UMA1000

Dimensions: 190x100x80mm  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Case Size: 40

R W S E HP LF



Slipper Pan -  
Female urinal

Code: SPN

Dimensions: 335x245mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Case Size: 25

R W S E HP LF



Urinal Bottle -  
Female urinal with handle

Code: UF1000

Dimensions: 170x100x80mm  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Case Size: 30

R W S E HP LF



Urinal Pan -  
Female Urinal

Code: PUF1000

Dimensions: 165x85x30mm  
Volume: 1000ml  
Material: Polyethylene  
Colour: White  
Case Size: 20

R W S E HP LF



Non Leak Valve -  
For use with UMA1000 &  
UML1000

Code: UMT

Dimensions: 50x90mm  
Volume: N/A  
Material: Polypropylene  
Colour: White  
Case Size: 25

R W S E HP LF



Urinal Bottle with UMT -Male urinal  
bottle with non leak valve  
Graduated

Code: USP1200

Dimensions: 180x100x85mm  
Volume: 1000ml  
Material: Polypropylene and  
Polyethylene  
Colour: N/A  
Case Size: 10

R W S E HP LF



Urinal Bottle Holder -  
For UBS1000, UML1000 and  
UM1000

Code: UBH2002

Dimensions: N/A  
Volume: 1200ml  
Material: Plastic Coated Wire  
Colour: White  
Case Size: 10

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

# Drinking Beakers, Cups, Lids and Denture Cups



Beaker - Feeder/drinking cup  
Graduated

Code: BFC250

Dimensions: 62x75mm  
Volume: 250ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 300



Beaker - Feeder/drinking cup with  
handles

Graduated  
Code: BFH250

Dimensions: 62x80mm  
Volume: 250ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 140



Lid -  
Narrow Spout for BFC250/  
BFH250

Code: BNS

Dimensions: 4x4mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Case Size: 300



Lid -  
Medium Spout for BFC250/  
BFH250

Code: BMS

Dimensions: 7x4mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Case Size: 300



Lid -  
Wide spout for BFC250  
& BFH250

Code: BWS

Dimensions: 13x8mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Case Size: 300



Lid -  
Rim feeder for beakers BFC250 &  
BFH250

Code: BRF

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Case Size: 200



Drinking Cup -  
With handles, spout lid and feeder  
lid

Code: GPMC6250

Dimensions: 73x50mm  
Volume: 250ml  
Material: Polypropylene  
Colour: Cream  
Case Size: 90



Teapot Style Beaker

Code: TFB

Dimensions: 88x55mm  
Volume: 250ml  
Material: Polypropylene  
Colour: White  
Case Size: 45



Beaker Set - Feeder/drinking cup  
with handle, spout lid and feeder lid  
Graduated

Code: HDMS

Dimensions: 74x90mm  
Volume: 300ml  
Material: SAN  
Colour: Clear  
Case Size: 40



Clear Glass Style Beaker

Code: TBC

Dimensions: 70x70mm  
Volume: 280ml  
Material: SAN  
Colour: Clear  
Case Size: 120



Denture Cup -  
With strainer and lid

Code: DNT

Dimensions: 90x80x60mm  
Volume: N/A  
Material: Polypropylene  
Colour: Green  
Case Size: 75



R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

# Jugs and Jug Sets



Jug  
Graduated

Code: JG500

Dimensions: N/A  
Volume: 500ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 50



Jug  
Graduated

Code: JG1000

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 45



Jug with lid  
Graduated

Code: JGL1000

Dimensions: 120x125mm  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 40



Jug Set - Includes tray, carafe with lid and beaker  
Graduated

Code: LST

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: N/A  
Case Size: 20



# Medicine Measures



Medicine Measure -  
Graduated

Code: MM60

Dimensions: 48x45mm  
Volume: 60ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 1500



Medicine Measure -  
Graduated black lines

Code: MM60BL

Dimensions: 48x45mm  
Volume: 60ml  
Material: Polypropylene  
Colour: Natural  
Case Size: 1500



Distributed By



# Patient Trays, Kidney Dishes and Vomit Bowls



Dispensing/Dressing Tray

**Code: ITC2715**

Dimensions: 270x150x30mm  
 Volume: N/A  
 Material: Polypropylene  
 Colour: Green  
 Case Size: 60

**R W S E HP LF**



Kidney Dish  
 Graduated

**Code: KDC250**

Dimensions: 250x120x50mm  
 Volume: 750ml  
 Material: Polypropylene  
 Colour: Green  
 Case Size: 60

**R W S E HP LF**



Vomit Bowl - With handle  
 Graduated

**Code: VB300**

Dimensions: 300x145x60mm  
 Volume: 1500ml  
 Material: Polypropylene  
 Colour: Green  
 Case Size: 25

**R W S E HP LF**

**R** Reusable

**W** Washer Decontamination

**S** Suitable for Steam Sterilization

**E** Suitable for Ethylene Oxide Sterilization

**HP** Suitable for Hydrogen Peroxide Sterilization

**LF** Latex Free

